

## OPTIONAL ACCESSORIES

### DENSITY BALANCES

0.1 mg - 10 kg

- Density Determination Kit is available in 2 types one to be mounted on the pan and the other type is with attachment for weigh below the balance. Density Determination Kit consists of pan (for weighing in air), suspended spiral pan (for weighing in water), S.S. stand for holding pan fitted to S.S. pan placed on the balance, stainless steel stand for placing beaker, glass beaker, optional sinker for testing materials having density less than water and software (program) for density determination.
- Facility for programming density of liquid (other than water) and temperature of water.
- Incase of multiple samples, storage of all the sample weights in air first and then in water, viewing of density of individual samples one by one on the display (Optional).
- Determination of density of liquid samples (Optional).
- Determination of air release value for oil samples (Optional).

