

ACCESSORIES

Magnetic Stirrer MS-01

Stirring speed : 0~2300 rpm at no load.
 Base Diameter : 120mm
 Maximum stirring capacity : 1000ml
 Power supply : 6V DC Adaptor
 Dimension : 120 dia x 50mm height
 Meter Weight : 500 gms.



Electrode Holder EH-01

Dimensions (mm) : 115(L) x 35(W) x 28(H)
 Weight : 22 g
 Suitable for Electrode : Dia 13mm~16mm



Stirr Bead SB-01

Dimensions : 27(L) x 10(W) x 6(H)
 Weight : 5g



Mini Printer

CMP-01

Specifications

CMP-02



Thermal
CMP-01

Thermal :	Printing Method	: Dot Matrix
Alpha Numeric :	Character Type	: 8x16dots
30 max. :	Character per Line	: 24
RS 232 :	Communication	: RS 232
9600 :	Baud Rate	: 9600
90mm per sec. :	Printing Speed	: 0.7 line per sec.
57mm :	Paper Width	: 57mm
8v,3A DC :	Power Supply	: 5v DC, 1.5 A
135(L) x 195(D) x 100(H) :	Dimension (mm)	: 110(W) x 180(H) x 90(D)
800g approx. :	Net Weight	: 600g approx.



Dot Matrix
CMP-02

Power Supply SMPS-04

Input : 100~300 V AC, 47-63 Hz
 Output : 9VDC, 500mA
 Dimensions with pin (mm) : 50(L) x 70(H) x 70(D)
 Dimensions without pin (mm) : 50(L) x 70(H) x 37(D)
 Weight : 100g



Battery Backup BBEC-01

Input : 100~300 V AC, 47-63 Hz
 Output : 9VDC, 500mA
 Battery : 6V, 1.2 Ah
 Deminsions : 150(L) x 65 (w) x 40(H)
 Weight : 600g



Calibration Solutions

Solution	Corresponding Electrode	Description	Accuracy
CSPH 4.0	pH	pH 4.0	± 0.1 pH
CSPH 7.0		pH 7.0	± 0.1 pH
CSPH 9.2		pH 9.2	± 0.1 pH
CSCN 1413	Conductivity	1413 µs/cm	±0.5%

