



BALANCE ACCESSORIES

BREEZE SHIELD

Contech provides glass breeze shield as standard with EMFC balances having accuracy of 1mg or better to shield the weighing area from breeze. For balances with accuracy less than 1mg, breeze shield can be provided as optional. Different types of Breeze Shields are available based on the model and application.



CONTECH[®]
Instruments Ltd.

AN ISO 9001, 14001 & 45001 COMPANY

Where innovation never stops.....

301, Punit Indl. Premises, Turbhe Naka, Navi Mumbai - 400 705.

Tel. : +91 22 6139 3000 (100 lines)

E-mail : sales@contechindia.in

Website: www.contechindia.com

FACTORY / NABL LAB



Plot No. EL-221 TTC Indl. Area, MIDC (Electronic Zone),
Mhape, Navi Mumbai - 400 710.