

BABY WEIGHING SCALES

5g to 20Kg


Features :

- ♦ Overload and shock load protection.
- ♦ Bright LED display.
- ♦ Large Pan Size.
- ♦ Pan Electrically Insulated from the scale for safety.
- ♦ Optional battery backup.
- ♦ Display hold facility.


Specifications :

Models	CT-20K5
Capacity	20kg
Readability	5g
Pan Size (mm)	490x280

Taring Range : Up to full capacity
 Response time : 3-5 secs
 Power Req : 90-250V, 50/60Hz, 10 Watt
 Operating Temperature : 10 to 45 deg C



TOP VIEW



FRONT VIEW

LEVELING FOOT
4-LEGS

All Specifications and designs are subject to change for improvement.

